# LaboForce





Specimen movers
for semi-automatic
preparation of
materialographic
specimens on LaboPol
grinding and polishing
machines

### The LaboForce specimen movers offer you:

- Semi-automatic fine grinding and polishing of single specimens
- Quick and easy operation
- High quality finished specimens
- Adjustable drip lubricator

The LaboForce specimen movers are perfect for small laboratories with modest specimen throughput and equally ideal for large labs, where the need to prepare single specimens occurs from time to time.

### Integration with LaboPol

LaboForce-3 and LaboForce-1 are prepared for use with the Struers grinding, lapping and polishing machines LaboPol-1, LaboPol-2, LaboPol-4, LaboPol-5, LaboPol-6 and LaboPol-21.

### Quick and easy operation

As both LaboPol and LaboForce are designed for maximum compatibility a very user-friendly operation is obtained. LaboForce is connected directly to the LaboPol, and all start and stop functions are controlled via the LaboPol, resulting in a fast and safe operation.

## Preparation of single specimens results in consumable and time savings

Both LaboForce-3 and LaboForce-1 are designed for the preparation of single specimens. The load on each specimen is applied individually. That makes

it possible to run 1, 2 or 3 specimens at a time without any need for blank or "dummy" specimens. It is also not necessary to plane grind the specimens as they are levelled individually. Both the preparation of single specimens and the exclusion of the plane grinding step results in considerable savings of consumables and preparation time.

### LaboForce-3

LaboForce-3 is running with a speed of 250 rpm and designed for fine grinding and polishing of 1-3 single specimens. The specimens are easily placed in the specimen mover plate and the force is set individually, adjustable from 5-40 N. Fine grinding can be carried out using MD-Allegro, MD-Largo or a single step on a fine grained SiC-Paper or diamond grinding disc. For polishing, Struers offers a wide range of cloths.

The specimen mover plates can be exchanged quickly to accommodate 25 mm (1"), 30 mm, 1¼", 1½" and 40 mm dia. specimens.

### LaboForce-1

LaboForce-1 is the economic choice. It is running with 8 rpm and suitable for fine grinding and polishing of 1-3 single specimens. Specimen sizes from 25 mm (1"), 30 mm, 1½", 1½", and 40 mm dia. are accommodated. The force can be set to between 2 - 20 N. Like with LaboForce-3, MD-Allegro, MD-Largo, fine grained SiC-paper or diamond grinding discs can be used for



fine grinding. The large range of Struers polishing cloths can be used.

### Adjustable drip lubricator

Part of the standard equipment for both LaboForce-3 and LaboForce-1 is an adjustable and exchangeable drip lubricator. This enables quick changes between different lubricants.

The dosing level can be accurately adjusted by a screw, ensuring reproducible preparation results and optimal use of lubricant.



LaboForce can also be used on the two disc grinding and polishing machine LaboPol-21



With LaboPol/LaboForce, the operator may lap, grind or polish as required by a quick and simple exchange of discs

Technical Data		
	LaboForce-3	LaboForce-1
Motor	1 Phase	1 Phase
Output	25 W	0.5 W
Amperage	220 - 240 V, 0.27 A	220 - 240 V, 16 mA
	100 - 120 V, 0.55 A	100 - 120 V, 32 mA
	Motor with thermal switch	
Rotational speed	250 rpm at 50/60 Hz	8/9.5 rpm at 50/60 Hz
Force per specimen	Min./max. approximately	Min./max. approximately
	5 N / 40 N	2 N / 20 N
Dimensions:		
Height	320 mm/12.6"	295 mm/116"
	(above grinding/polishing disc)	(above grinding/polishing disc)
Width	185 mm/7.5"	185 mm/7.5"
	(including drip lubricator)	(including drip lubricator)
Depth	325 mm/12.8"	265 mm/10.5"
Weight	9.5 kg/21 lbs	5.9 kg/13 lbs
Safety standards	General: 89/392/EEC & 91/386/EEC & 93/44/EEC EN292.1, EN292-2, ENB60204-1 (IEC204-1) EMC: 89/336/EEC & 92/231EECEN50081-1, EN50082-1 FCC part 15, subpart B	

Safety standards	General: 89/392/EEC & 91/386/EEC & 93/44/EEC EN292.1, EN292-2, ENB60204-1 (IEC204-1) EMC: 89/336/EEC & 92/231EECEN50081-1, EN50082-1 FCC part 15, subpart B		
Specifications		Code	
operated from LaboPol-	en mover for the preparation of 1-3 specimens. Is mounted on and -21, LaboPol-6, LaboPol-5, LaboPol-2 or LaboPol-1. Complete with cimen mover plates and additional drip lubricators are ordered	LAFO3	
Specimen mover plate for 3 x Ø 25 mm or 1" s		LAFTY	
for 3 x Ø 30 mm specimens		LAFET	
for 3 x Ø 1 1/4" specimens		LAFQU	
for 3 x Ø 1 ½" specimens		LAFHA	
for 3 x Ø 40 mm specimens		LAFYF	
and operated from Labo	en mover for the preparation of 1 - 3 specimens. Is mounted on pPol-21, LaboPol-5, LaboPol-4, LaboPol-2 or LaboPol-1. Complete Specimen mover plates and additional drip lubricators are ordered	LAFO1	
Specimen mover plate		PEMFI	
for 3 x Ø 25 mm or 1" specimens for 3 x Ø 30 mm specimens		PEMTE	
for 3 x Ø 1 1/4" specimens		PEMQU	
for 3 x Ø 1 ½" specimens		PEMNA	
for 3 x Ø 40 mm specimens		PEMYF	
· ·	orForce-3 and LaboForce-1  or for quick changing between different lubricant types		

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice



Struers A/S Valhøjs Allé 176 DK-2610 Rødovre Phone +45 36 70 35 00 Fax +45 38 27 27 01 e-mail: struers@struers.dk www.struers.com



### DEUTSCHLAND

Struers GmbH Karl-Arnold-Strasse 13 B D - 47877 Willich Telefon +49/(0) 21 54/4 86-0 Telefax +49/(0) 21 54/4 86-222 e-mail: verkauf.struers@struers.de

### ÖSTERREICH

Struers GmbH

Zweigniederlassung Österreich Zweigniederiassung Osteri Ginzkeyplatz 10 5020 Salzburg Telefon (0662) 625711 Telefax (0662) 625711-78 e-mail: stefan.lintschinger@struers.de

### SCHWEIZ

**Struers GmbH** Zweigniederlassung Schweiz Weissenbrunnenstrasse 41 CH-8903 Birmensdorf Telefon (01) 77763-07 Telefax (01) 77763-09 e-mail: rudolf.weber@struers.de

### HOLLAND

Struers GmbH Nederland Electraweg 5 NL-3144 CB Maassluis Tel.: +31 (0) 10 599 72 09 Fax: +31 (0) 10 599 72 01 e-mail: glen.van.vugt@struers.de

### FRANCE

### Struers S.A.S.

370, rue du Marché Rollay F- 94507 Champigny sur Marne Cedex Téléphone +33 1 55 09 14 30/31 Télécopie +33 1 55 09 14 49 e-mail: struers@struers.fr

### BELGIQUE Struers S.A.S. 370, rue du Marché Rollay

F- 94507 Champigny sur Marne Cedex Téléphone +32 43 70 93 18 Télécopie +32 43 70 93 19 e-mail: struers@struers.fr

### JAPAN

### Marumoto Struers K.K.

Takara 3rd Building 18-6, Higashi Ueno 1-chome Taito-ku, Tokyo 110-0015, Phone: +81 3 5688-2914 Fax: +81 3 5688-2927 e-mail: struers@struers.co.jp

### SINGAPORE Struers A/S

10 Eunos Road 8 #12-06 North Lobby Singapore Post Centre Singapore 408600 Phone +65 6299 2268 Fax +65 6299 2661 e-mail: struers.sg@struers.dk

### UNITED KINGDOM Struers Ltd. Erskine Ferry Road

Old Kilpatrick Glasgow, G60 5EU Phone 01 389 877 222 Fax 01 389 877 600 e-mail: info@struers.co.uk

### USA

### Struers Inc. 810 Sharon Drive

Westlake, OH 44145 Phone (888) 787-8377 Fax (440) 871 8188 e-mail: info@struers.com